

Product Code . JA-BBE-3557



Blood Bank Centrifuge Counter Balance

Description

Blood Bank Centrifuge Counter Balance

Whole blood is loading of blood bags in high speed centrifuge machine.

Before loading the blood bag, the weights of the two loads (each blood bag + bucket) have to be balanced in such a way that their weights should be the same or very close to each other.

This is a vital step of plasma separation process to save expensive centrifuge machine.

By conventional methods the process of balancing is time consuming, inaccurate and cumbersome.

This equipment makes the whole process of balancing execute shortly and accurately.

It is the only way to save expensive centrifuge machine from damage.

Jain Laboratory Instruments Pvt. Ltd,
Hargolal Road, Ambala Cantt, Haryana India
Direct Contact Details ☎ +91-8569909696 ✉ sales@jlabexport.com
💻 www.jlabexport.com