

JLab Instruments







Product Code . JA-O-3899

Forced Convection Laboratory Ovens

Description

Forced Convection Laboratory Ovens

Safe, Superior Protection for Sample, User and the Environment.

Multiple redundant over-temperature protection systems to guarantee maximum sample and user protection.

Electronic over-temperature protection built into the microprocessor.

Redundant mechanical over-temperature protection, adjustable from the front, independent from the microprocessor.

Overall temperature protection meets.

Red LED illuminates if external mechanical temperature protection is engaged.

Controller will control temperature at the over temperature set point.

All electrical components UL recognized.

Electrical circuit protection in accordance with UL requirements.

